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---

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1. Chen CH, Nevo E, Fetics B, *et al.* Estimation of central aortic pressure waveform by mathematical transformation of radial artery tonometry: validation of generalised transfer function. *Circulation* 1997;95:1827-36.
2. Karamanoglou M, O'Rourke MF, Avolio AP, *et al.* An analysis of the relationship between central aortic and peripheral upper limb pressure waves in man. *Eur Heart J* 1993;14:160-7.
3. Pauca AL, O'Rourke MF, Kon ND. Prospective evaluation of a method for estimating ascending aortic pressure from the radial artery pressure waveform. *Hypertension* 2001;38:932-7.
4. Wilkinson IB, Fuchs SA, Jansen, IM, *et al.* Reproducibility of pulse wave velocity and augmentation index measured by pulse wave analysis. *J Hypertension* 1998;16:2079-84.
5. Seibenhofner A, Kemp CRW, Sutton AJ, *et al.* The reproducibility of central aortic blood pressure measurements in healthy subjects using applanation tonometry and sphygmocardiography. *J Hum Hypertens* 1999;13:625-9.
6. Filipovsky J, Svobodova V, Pecan L. Reproducibility of radial pulse wave analysis in healthy subjects. *J Hypertension* 2000;18:1033-40.
7. Savage MT, Ferro CJ, Pinder SJ, *et al.* Reproducibility of derived central arterial waveforms in patients with chronic renal failure. *Clin Sci* 2002;103:59-65.