



# The Role of Central Aortic Pressure in the Management of High Blood Pressure

**CASE STUDIES**

# Objective of Presentation

**The presentation will outline hypothetical cases that illustrate the potential utility of PWA for management decisions based on the following:**

- Clinical evidence supports incorporating non-invasive central aortic blood pressure measurements into blood pressure management.
- Numeric thresholds for central BP variables can be defined utilizing published research.
- The SphygmoCor<sup>®</sup> System can help guide treatment decisions designed to prevent or reduce long-term target organ damage and cardiovascular events resulting from increased aortic pressure and therefore reduce the socioeconomic burden of hypertension.

# Rationale for Utilizing PWA and Central BP Variables: Optimizing Hypertension Pharmacotherapy

## Pharmacotherapy Concerns

- Undertreatment, overtreatment, compliance, drug cost, adverse events, drug interactions
- Generally, lifetime treatment
- All above impact a patient's adherence behavior to prescribed treatment and the burden of hypertension

## Prescription Optimization

- Optimizing prescription medication including self-administration of therapy is critical to controlling hypertension

# Rationale for Utilizing PWA and Central BP Variables: Optimizing Hypertension Pharmacotherapy

- Incorporation of non-invasive measurements of central aortic pressures can improve hypertension management in the following areas:
- Refine monitoring requirements.
- Reduce over-treatment.
- Improve under-treatment.
- Reduce costs of management (e.g., medication costs, monitoring such as ambulatory blood pressure monitoring (ABPM)).

# Rationale for Utilizing PWA and Central BP Variables: Optimizing Hypertension Pharmacotherapy

1. Incorporation of PWA into the treatment paradigm for hypertension (i.e., in addition to brachial pressure monitoring) can provide clinically relevant value to patient care:
2. Confirmation of hypertension so that initiation of medication is more likely to be the correct decision for an individual patient  
**Scenario:** *Concurrent elevation in brachial and central pressures*
3. Avoiding initiation of medication when white coat hypertension is suspected  
**Scenario:** *Elevated brachial pressure and normal central pressures*
4. Confirmation that increased treatment may not be needed  
**Scenario:** *Borderline high peripheral pressures and normal central pressures*
5. Targeting when to consider reduction of medication  
**Scenario:** *Normal peripheral and low central pressures, or extended period of normal peripheral and normal central pressures (particularly in the setting of medication tolerance issues)*

# Central Aortic Pressures: Clinical Thresholds\*

## Central aortic systolic pressure

- Normal: < 112 mmHg
- Elevated: 113 – 130 mmHg (equivalent to Stage 1 hypertension)
- High: >130 mmHg (equivalent to Stage 2 hypertension)

## Augmentation Index (Aix)

- Normal: < 35%
- High: > 35%

## Central aortic pulse pressure

- Normal: <50 mmHg
- High: > 50 mmHg

## Central aortic diastolic pressure

- Expected to be similar to brachial diastolic pressure. Threshold for brachial pressure can be used.

\*based on research and review peer-reviewed publications.

# Case Examples\*

\*Cases are not actual patients and have been developed to illustrate utility of PWA

# Clinical Cases

**Case 1:** Untreated hypertension

**Case 2:** Persistent systolic hypertension on drug therapy

**Case 3:** Suspicion of white coat hypertension

**Case 4:** Possible adverse effects of hypertension medications

**Case 5:** Reducing hypertension medications

# Case 1a: Untreated Hypertension

**History:** 48-years, male, moderately overweight, sedentary, non-smoker, no previous history of hypertension

**Reason for Visit:** annual routine assessment

## **Clinical Findings:**

- HR: 76 bpm, BMI: 29 kg/m<sup>2</sup>
- Brachial BP: 145/90 mmHg
- Central BP: 129/88 mmHg, Alx = 38%
- No evidence of secondary hypertension

## **Diagnosis:**

- Stage 2 hypertension, minimal elevation in central BP

## **Interpretation:**

- Brachial BP elevation either essential hypertension or white coat hypertension
- Given minor elevation in central BP, extended trial of lifestyle modification may be indicated prior to medication initiation.
- Follow-up: clinic visit for reassessment. Consider 24-hour ABPM.

# Case 1b: Untreated Hypertension

**History:** 56-years, female, physically active, current smoker 20 pack years, no previous history of hypertension

**Reason for Visit:** routine assessment

## **Clinical Findings:**

- HR: 78 bpm, BMI: 24 kg/m<sup>2</sup>
- Brachial BP: 148/92 mmHg
- Central BP: 138/90 mmHg, Alx = 48%
- No evidence of secondary hypertension

## **Diagnosis:**

- Stage 2 hypertension, elevation in central BP and Alx

## **Interpretation:**

- Elevation in both brachial BP, central BP and Alx. Increased risk for vascular end-organ damage.
- Prompt follow-up is indicated to confirm elevated BP. Confirmation should result in initiation of medication. Counseling for smoking cessation.
- Follow-up: Clinic visit for reassessment, recommend home BP monitoring.

# Case 2a: Persistent Systolic Hypertension on Drug Therapy

**History:** 62-years, female, physically active, ex-smoker (30 pack years), known essential hypertension, asymptomatic

**Reason for Visit:** BP follow-up

## **Clinical Findings:**

- HR: 64 bpm, BMI: 26 kg/m<sup>2</sup>
- Brachial BP: 130/86 mmHg
- Central BP: 116/84 mmHg, Alx = 15%

## **Diagnosis:**

- Stage 1 hypertension (on treatment), almost normal central variables

## **Interpretation:**

- Acceptable control of brachial BP although not normal. No indication of increased risk based on central BP and Alx.
- Although research suggests that more aggressive brachial BP control may be beneficial, the near normal central BP informs management decision towards maintaining current therapeutic measurement.
- Follow-up: Routine BP follow-up clinic visit as per usual schedule.

# Case 2b: Persistent Systolic Hypertension on Drug Therapy

**History:** 62-years, female, physically active, ex-smoker (30 pack years), known essential hypertension, asymptomatic

**Reason for Visit:** BP follow-up

## **Clinical Findings:**

- HR: 72 bpm, BMI: 26 kg/m<sup>2</sup>
- Brachial BP: 132/88 mmHg
- Central BP: 126/88 mmHg, Aix = 48%
- No evidence of secondary hypertension

## **Diagnosis:**

- Stage 1 hypertension (on treatment), elevated central variables

## **Interpretation:**

- Borderline acceptable control of brachial BP although above normal. Additional increased risk of adverse vascular events based on central BP and Aix.
- The combination of borderline acceptable brachial BP and elevated central BP variables suggests that more aggressive hypertension management is indicated, through increasing dose of current medications (if possible) and/or adding another class of medication.
- Follow-up: More frequent monitoring including home BP monitoring.

# Case 3: Suspicion of White Coat Hypertension

**History:** 32-years, male, mildly overweight, sedentary, non-smoker, hypercholesterolemia on routine testing

**Reason for Visit:** annual assessment

## **Clinical Findings:**

- HR: 82 bpm, BMI: 27 kg/m<sup>2</sup>
- Brachial BP: 145/86 mmHg
- Central BP: 112/84 mmHg, Alx = 20%

## **Diagnosis:**

- Stage 2 hypertension based on brachial BP, normal central BP and Alx

## **Interpretation:**

- Brachial BP elevation either essential hypertension or white coat hypertension (WCH)
- Suspect white coat hypertension based on normal central pressures. However, data suggest that white coat hypertension may increase risk for vascular adverse events that is between normal (no risk) and essential hypertension. Lifestyle counseling (exercise, diet, weight loss) and investigation for WCH is indicated.
- Follow-up: Clinic visit for reassessment. 24-hour ABPM.

# Case 4: Possible Adverse Effects of Hypertension Medications

**History:** 74-years, female, never-smoker, sedentary, diabetes (controlled), NYHA class II, LVEF 45%

**Reason for Visit:** increased dyspnea, fatigue, presyncope on standing

Medications: furosemide, oral hypoglycemics, ARB, CCB

## **Clinical Findings:**

- HR: 74 bpm, BMI: 32 kg/m<sup>2</sup>
- Brachial BP: 124/82 mmHg (sitting), 118/78 mmHg (standing)
- Central BP: 108/78 mmHg, AIx = 39%

## **Diagnosis:**

- Normal brachial hypertension, low central BP, mild elevation in AIx

## **Interpretation:**

- Controlled brachial BP. Low central BP with slight elevation in AIx.
- Low central BP suggests that hypotension may be contributing to symptoms.
- Consider reducing dose or stopping at least one medication for hypertension. Central BP provides additional data supporting safety of stopping one medication (e.g., central BP is not elevated and is low).
- Follow-up: Prompt clinic visits for reassessment and monitoring of BP and symptoms. Should be recommended for home BP monitoring if not currently using HBPM. Consider ABPM.

# Case 5a: Controlled Hypertension - Reducing Hypertension Medications

**History:** 56-years, female, smoker (40 pack-years), moderately active, diabetic (controlled on oral hypoglycemic drug)

**Reason for Visit:** annual assessment, treated for hypertension (diuretic, ACEI)

## **Clinical Findings:**

- HR: 68 bpm, BMI: 24 kg/m<sup>2</sup>
- Brachial BP: 124/84 mmHg
- Central BP: 118/80 mmHg, Alx = 48%

## **Diagnosis:**

- Brachial BP normal, mild elevation in central BP, high Alx

## **Interpretation:**

- High Alx indicates increased risk for vascular end-organ damage.
- Counseling for smoking cessation. Given high Alx, should not decrease medication but may consider replacement of a medication class to attempt to lower Alx.
- Follow-up: Clinic visit for reassessment, recommend HBPM.

# Case 5b: Controlled Hypertension - Reducing Hypertension Medications

**History:** 56-years, female, physically active, smoker 22 pack years

**Reason for Visit:** annual assessment, treated for hypertension (2 classes)

## **Clinical Findings:**

- HR: 75 bpm, BMI: 24 kg/m<sup>2</sup>
- Brachial BP: 124/84 mmHg
- Central BP: 118/72 mmHg, Alx = 18%

## **Diagnosis:**

- Brachial BP normal, mild elevation in central BP, normal Alx


## **Interpretation:**


- BP controlled. Can consider reducing medications to minimum needed to control BP
- Follow-up: Clinic visit for reassessment, recommend HBPM.

# How to Use Central Blood Aortic Blood Pressure

# Guidelines for Hypertension (2017)

2017 ACC/AHA/AAPA/ABC/ACPM/AGS/APhA/ASH/ASPC/NMA/PCNA Guideline for the Prevention, Detection, Evaluation, and Management of High Blood Pressure in Adults  
 A Report of the American College of Cardiology/American Heart Association Task Force on Clinical Practice Guidelines






## PROPOSED GUIDELINE

### Scenario: Normal Brachial BP

- No change in guideline


\*cSBP cannot exceed brachial systolic BP (i.e., will always be < 120 mm Hg)



## PROPOSED GUIDELINE

### Scenario: Elevated Brachial BP

- Diagnosis of White Coat Hypertension when cSBP is normal
- Schedule routine follow-up



## PROPOSED GUIDELINE

### Scenario: Stage 1 Hypertension

- Largest impact in management decisions
- Diagnosis of White Coat Hypertension when cSBP is normal
- cSBP hypertension – confirmation of hypertension and initiation of treatment
- Incorporate cSBP in follow-up for management decisions related to increasing or decreasing medications

# PROPOSED GUIDELINE

## Scenario: Stage 2 Hypertension

- Diagnosis of White Coat Hypertension when cSBP is normal
- Schedule follow-up in 3 to 6 months
- Incorporate cSBP in follow-up for management decisions related to increasing or decreasing medications

